


NEW PRODUCTS

**WATER AND SALINE NOW AVAILABLE
FROM CHRONO-LOG**



**INCORRECT SALINE
(Injection/Infusion)
... Inhibits Aggregation**

**INCORRECT WATER
(Millipore)
... Causes Abnormal
Luminescence**

**Chrono-log Water and Saline Do Not Contain Preservatives/
Additives and provide:**

- **Improved Accuracy**
- **Reduced Costs**
- **Time, Sample and Reagent Savings**

-> P/N 397 - Saline 0.9% (15 mL)

-> P/N 398 - Purified Water (15 mL)

**TO IMPROVE YOUR PLATELET TESTING, ORDER
PRESERVATIVE/ADDITIVE-FREE WATER AND SALINE FROM:**

CHRONO-LOG CORP.

2 WEST PARK ROAD, HAVERTOWN, PA 19083-4691, USA

Phone: 1-800-247-6665; In PA 610-853-1130; Fax 610-853-3972

email: chronolog@chronolog.com

www.chronolog.com

